

Microfluidic chip with QCM

for detection of influenza virus subtype

OS. Kagiyama¹ T. Masuda¹ K. Kaihatsu² N. Kato² F. Arai¹



¹Nagoya University, JAPAN

²Osaka University, JAPAN



NAGOYA UNIVERSITY

New method for detection of influenza subtype

Background

Infectious disease


- Influenza
- Dengue fever etc...



Need simple and accurate Diagnosis

Previous method

- Rapid diagnosis kit
- **Rapidly** **Specifically**




Conventional technique

- Virus purification by Hydroxyapatite
- Viral DNA/RNA extraction by Silica


Purpose

Establish the method to detect viral RNA.


Concept





Peptide nucleic acid




Sensitizer
(Gold colloids)



1. Viral RNA captured **specifically** by Peptide Nucleic Acid (PNA).
2. Gold colloids captured to **increase sensing mass**.



Sensitizer




	Mass [g/piece]	Size [nm]
Viral RNA	1.2×10^{-18}	40~80
Gold colloid	6.5×10^{-16}	$\Phi 40$


540 times heavier

Experiment

Setting




Detection of viral RNA




Design

QCM




Thickness : 0.1 mm

Microfluidic channel



Thickness : 0.1 mm



Reference

- [1] M.Niimi et al, "Virus purification and enrichment by hydroxyapatite chromatography on a chip", volume 201, 2014, P185-190 Sensors and Actuators B
- [2] K.Kaihatsu et al, "Sequence-Specific and Visual Identification of the Influenza Virus NS Gene by Azobenzene-Tethered Bis-Peptide Nucleic Acid" 2013, PLOS ONE
- [3] S.Kagiyama et al, "Detection of influenza virus subtype using quartz crystal microbalance" 2014, Micro-Nano Mechatronics and Human Science (MHS)

Contact person: Taisuke Masuda

E-mail: masuda@mech.nagoya-u.ac.jp, URL: http://www.biorobotics.mech.nagoya-u.ac.jp/

Furo-cho, Chikusa-ku, Nagoya 464-8603, Japan

Dept. of Micro-Nano Systems Engineering, Nagoya University

TEL: 052-789-5026, FAX : 052-789-5027

Acknowledgements :

This work was supported in part by the Nano-Macro Materials, Devices and System Research Alliance Project from the Ministry of Education, Culture, Sports, Science and Technology of Japan (MEXT).

	Rapid diagnosis kit (BD Directigen TM FluA)	Our chip
Detection limit	$10^3 \sim 10^4$	5×10^3
Detect subtype	x	○